



# Maintain Peak Performance with Wearable Laser Therapy

Get on board with this convenient, portable, wearable and safe method of Low Level Light Therapy (LLLT). Spectra Therapy FLASHwrap™

- Reduces pain & inflammation
- Improves circulation
- Offers quicker healing with less infection
- Apply to any area of your horse

**Wrap dispenses treatments  
AUTOMATICALLY.  
Programmable from every  
10 minutes - 24 hours.**



## PAIN RELIEF FOR YOUR HORSE

Laser therapy is recommended and used by many veterinarians. Spectra Therapy's FLASHwrap™ is one of the most innovative and promising laser therapy devices released to date. It is aimed and approved to provide micro stimulation of cells indicated for the temporary relief of muscle and joint aches and pains, arthritis and muscle spasms, relieving stiffness, promoting relaxation of muscle tissue and to temporarily increase blood circulation.

Within minutes, this automated unit delivers specific rates of therapy frequencies and will time out as programmed which can be prescribed by a veterinarian. FLASHwrap™ laser therapy is most effective in reducing edema. LED products that generate heat are the worst thing for swollen areas before working out. Laser therapy is the most effective.

Unlike cumbersome units which require one to hold the device and move it manually around the area for 10 to 20 minutes, this unit allows one to wrap it around any affected area providing bursts of therapeutic laser rays of healing light into cells, tissues and arteries of your horse.







# Reduce Inflammation with Wearable Laser Technology

## HOW IT WORKS:

Light energy is emitted from the laser in the form of photons into the damaged cells. As the photons are absorbed by the cells through laser therapy it stimulates the mitochondria to accelerate production of ATP [Adenosine Triphosphate].

This biochemical increase in cell energy is used to help transform cells from a state of illness to a stable healthy state.

## BIOLOGICAL BENEFITS OF THERAPEUTIC LASERS

Laser therapy is used by millions all over the world to live a better life with less pain and better blood circulation. It is a noninvasive, drug free, safe means to deliver photon energy to specific treatment areas and depths without risk or harm to you, the patient.

Anti-inflammatory action laser light reduces swelling caused by bruising or inflammation of joints to give enhanced joint mobility.

Faster wound healing fibroblast development and collagen synthesis in damaged tissue is accelerated by laser light stimulation.

Increased metabolic activity increases output of specific enzymes to provide greater oxygenation of blood cells and improved production of the basic food source for cells generating Adenosine Triphosphate (ATP).

Increased vascular activity laser light induces temporary vasodilation, increasing blood flow to damaged areas.

Rapid cell growth and cellular reproduction is accelerated by laser therapy. Reduction of fibrous tissue decreases scar tissue damage resulting from cuts, scratches, burns or post-surgical wounds.

## EASE OF USE

Spectra Therapy's equine-specific FLASHwrap™ is convenient to use, takes only a few minutes to apply and covers a larger area than comparable hand-held devices. The therapy can be applied to virtually any area of the horse -- joints, back, tendons, thighs, hocks, neck and hoof. Laser therapy is transmitted through filaments into the tissue.

Wrap before the performance event, delivering a few treatments to reduce joint pain and improve the oxygenation of the blood supply pre-charging the muscular and skeletal systems. Then wrap after the event to assist the lymph and adrenal systems for quicker healing with less infection and reduced pain. It is a one size fits all device, so you can use it on any of your horses.

***FLASHwrap™ should be a staple in every stable.***



# Be A Winner!



**Discover the healing power  
of equine-specific wearable laser therapy!**

- Easy • Safe • Portable • Effective

Call today to order your  
Spectra Therapy FLASHWRAP™ for your stable.

MADE IN THE USA

248 524-6300

Products Shipped Internationally • Distributor inquiries welcomed  
Engineering & manufacturing located in Troy & Shelby, Michigan

DISTRIBUTED BY:

